



Baycox Pumps, Denmark

Article No. and Description	Picture
<p style="text-align: center;">800148</p> <p>Baycox vet. pump for piglets only. Not adjustable (0,7 ml). Can be used on 250 and 1000 ml bottles but is mainly used with 250 ml bottles</p>	
<p style="text-align: center;">207265</p> <p>Baycox Bovis adjustable pump with hose. To be used in cattle. Adjustable between 3 and 30 ml. Can be used on 250 ml and 1000 ml bottles</p>	
<p style="text-align: center;">212593</p> <p>Baycox vet./Sheep adjustable pump without hose. Bottle is mounted on the top of the pump. To be used with piglets or lambs. Adjustable within the range of 0,2 – 2,0 ml. Can only be used with 250 ml bottles.</p>	